





READING THE ROBERTO REPORTS

If you're seeing scores in green with an occasional red or blue, that's normal. If scores are frequently red or blue, that means your performance is improving.




If you notice scores in red occurring frequently, it could signal there is an issue that should be addressed. Because each module targets a specific aspect of neurocognitive and neuromotor performance, data generated by Roberto can provide additional data for caregivers, parents, coaches, technicians and more.



FIND OUT MORE AT WWW.ROBERTOAPP.COM

TREND REPORTS AT YOUR FINGERTIPS

Roberto collects your past sessions for easy review. With continued use, you can log your neurocognitive and neuromotor health data. With that information, you can be proactive and make more informed decisions about your brain health.



Questions about the Roberto App or how to interpret the reports?
Email support@rc21x.com or
call 855-697-2219 x114

**The 6-Minute
Brain Performance
Program**

Track your
performance
on any
mobile device,
available
anywhere 24/7

RC21X
425 Mill Street, Suite 1
Coraopolis, PA 15108
WWW.RC21X.COM
WWW.ROBERTOAPP.COM



A STEP-BY-STEP GUIDE TO USING ROBERTO

The Roberto App uses fun, game-like experiences to monitor and track your brain performance. Here's how to use this innovative app and to begin gathering data that can be shared with family, coaches, trainers, physicians and more.

1

DOWNLOAD THE ROBERTO APP


Roberto works on any iOS or Android device and can be found in your device's app store.

2

REGISTER YOUR ACCOUNT

Each user must have a unique email address associated with his or her account. You'll be asked a series of questions, including name and date of birth. You can also choose to sign in with your Facebook account. **Please use the code WCFF2020.**

Be sure to enter the code **WCFF2020** when registering your account.



3

TEST YOUR DEVICE AND WATCH THE INSTRUCTIONAL VIDEO


Make sure your volume is turned on and that you can hear the test message. The video gives a quick overview of the app.

DOWNLOAD THE APP NOW
WWW.ROBERTOAPP.COM

4

START USING ROBERTO

Find a quiet space and allot at least 6 minutes of uninterrupted time. Follow the instructions on the screens before each module. Remember, Roberto rewards speed and accuracy. You'll receive your first Brain Performance Profile report within seconds after you're finished. After a few sessions, you have the option to skip the instructions.




In this module, we flash five shapes which must be repeated in the order seen.

5

ESTABLISH YOUR NORMAL RANGES


After four sessions, Roberto gets to know you. We call this earning your colors as all future Brain Performance Profile reports will be color-coded. Scores in green mean you're performing at your normal range. Scores in blue indicate a new high score. Scores in red mean you are below what's normal for you.



By "earning your colors," you can track your performance at a glance with our easy-to-read report.

ESTABLISH YOUR **NORMAL RANGE**.
MONITOR YOUR SCORES.
TRACK **LIFESTYLE CHANGES**.

In-app messaging reminds you to take your Roberto sessions regularly.



6

MAKE ROBERTO PART OF YOUR ROUTINE


We designed Roberto sessions to be quick for daily use. For kids playing sports, it's important to take one session after a game and another during the week.

By establishing **trend data** on your **brain performance**, you have **objective data** to share with a parent, caregiver or clinician in the event of **illness, injury or recovery**.

7

START DEVELOPING TREND DATA

Roberto records overall scores and scores within each module. You can also add notes in the report section. Those notes allow you to collect more detailed information to associate with your session. Maybe you slept poorly or were under stress when you used the app.



The notes feature allows you to add context by creating session comments for later review.